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Cambridge International Stem Cell Symposium

19th - 21st September 2018


206 The "Helene Medium": specialized stem cell culture medium

Yang CF, Saito M, Shirakawa K, Matsuo T


STEMCELL Co Ltd, Japan

Keywords: cell culture, transplantation


The stem cell culture requires many different conditions from normal or cancer cell culture, such as serum concentration, culture methods, and most of all, the culture medium. In our clinic in Japan, we provide stem cell treatments to our patients. We need more efficient and safer stem cell culture; thus, we have developed a specialized stem cell culture medium for primary stem cell culture. The "Helene Medium" is designed specially for stem cell growth, with better growth rate, stable cell growth, less chance of differentiation. We grow the stem cells in our medium and two other commercial mediums, examine the growth rate, cell morphology and passage number to evaluate the cell quality. Besides, serum usage in stem cell culture might lead to cell differentiation. We also test the different concentrations and types of serum. Stem cells could grow easily in Helene Medium but other cell types such as fibroblast cells, are difficult to grow. In our research, we have found that compared with other commercial mediums, stem cells cultured by our medium grow faster and with more stable cell morphology. Also, the expression of clusters of differentiation (CD) shows that stem cells are able to keep their potency during cell culture.


No significant difference is observed in either cases on Day 1. Whereas on Day 3, Helene Medium outstands the commercial mediums in terms of better cell growth, more stable morphology under both serum-free and serum-containing condition.

Stem cells cultured without serum grows slower than those with serum.

Yet in either cases, number of cells cultured in Helene medium exceeded that in the two commercial medium, showing that cells in Helene exhibits better growth rate.

Size of stem cells is considered as one of the validations for differentiation tendency, as stem cells are smaller than its differentiation product, fibroblasts (10-15 µm). From the result (left), only stem cells cultured within Helene medium demonstrates stem cells' character of smallness.

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